



DTR Medical

your
single-use
solution

General & Vascular Sterile Single-use Products

GENERAL & VASCULAR



www.dtrmedical.com



DTR Medical

General Benefits

- *Sterile for extra confidence*
- *High quality, robust & well balanced instruments suitable for all*
- *Cost effective alternative to re-processing*
- *Removes costly repairs or replacements*
- *Suitable for high risk procedures*
- *Supplied with peel off labels for full traceability*
- *Sterilised by gamma irradiation to minimise the risk of moisture or residues*
- *The range is developed with clinicians to replicate or better their re-usable instruments.*

General & Vascular Single-use Products

2011 Edition

Contents

Forceps other and Scissors	4
Forceps Bipolar	7
General and Vascular	8
Probes	11
Suction Handles	12
Tip Shields	16
Tubing Products	16

DTR Medical Ltd

17 Clarion Court, Enterprise Park, Swansea, SA6 8RF, UK
Telephone: +44 (0)1792 797910 Fax: +44 (0)1792 797955
Email: customerservice@dtrmedical.com


All products are sterile single-use unless otherwise stated

Forceps Other & Scissors

- Selected high quality single-use instruments
- Supplied sterile with traceability labels
- Compliments other single-use items


● CRILE & MOSQUITO

- CFS2001-D Crile Artery Forceps, 140mm Straight, double packed 20 units
- CFC2001-D Crile Artery Forceps, 140mm Curved on Flat, double packed 20 units
- MQT5001-S Mosquito Forceps, single packed 50 units
- MQC2012 Micro Mosquito Forceps 120mm Curved, double packed 20 units




● SPENCER WELLS

- SWF2010 Spencer Wells Artery Forceps 120mm straight, double packed 20 units
- SWF2015 Spencer Wells Artery Forceps 150mm straight, double packed 20 units
- SWF2018 Spencer Wells Artery Forceps 180mm straight, double packed 20 units
- SWF2020 Spencer Wells Artery Forceps 200mm straight, double packed 20 units




3 TUBING CLAMPS

- TCF2051 Tubing Clamps 150mm, double packed 20 units



ADSON & McINDOE

- ADS5001-S Adson Non Toothed Forceps, single packed 50 units
- ATS5001-S Adson Toothed Forceps, single packed 50 units
- MDF5001-S McIndoe Dissecting Forceps, 150mm non toothed, single packed 50 units



FINE SUTURE

- SPF2001-S Fine Suture Forceps, Plain 105mm, single packed 20 units
- SPF2001-D Fine Suture Forceps, Plain 105mm, double packed 20 units
- SSF2001-S Fine Suture Forceps, Serrated 105mm, single packed 20 units
- SSF2001-D Fine Suture Forceps, Serrated 105mm, double packed 20 units



RAMPLEYS


- RSH2181 Rampley Sponge Holders 180mm, double packed 20 units
- RSH2241 Rampley Sponge Holders 240mm, double packed 20 units



NEW

NEEDLE HOLDERS

- MNH2013 Mayo-Hegar Needle Holder 130mm, double packed 20 units
- MNH2015 Mayo-Hegar Needle Holder 150mm, double packed 20 units
- KNH2013 Kilner Needle Holder 135mm, double packed 20 units



NEW Range

● SCISSORS

- Stainless steel sharp sharp scissors
- First time sharpness which cut even at the tips
- Particularly suitable for clinic use
- Straight and curved options

Kilner

KSS2001-S Kilner Scissors, single packed 20 units

Iris

ISS2111-S Iris Scissors, Straight, Sharp Point, single packed 20 units

ISC2111-S Iris Scissors, Curved, Sharp Point, single packed 20 units

FSP2002-S Fine Suture Pack, single packed 20 units.

Stevens Tenotomy

TSS2001 Stevens Tenotomy Scissors, 115mm, single packed 20 units

Mayo

MSS2141 Mayo Scissor, Straight, Bevelled Blade, 140mm, double packed 20 units


MSC2141 Mayo Scissor, Curved, Bevelled Blade, 140mm, double packed 20 units

MSS2171 Mayo Scissor, Straight, Bevelled Blade, 170mm, double packed 20 units


MSC2171 Mayo Scissor, Curved, Bevelled Blade, 170mm, double packed 20 units



Kilner



Iris



Stevens Tenotomy



Mayo

NEW Range

SUTURE PACKS

- Consistently high quality instruments
- Reliable service
- Customised options available

- SPIA2003-S Suture Pack inc Iris Scissors, Adson Forceps and Mosquito, single packed 20 units
- SPKA2002-S Suture Pack inc Kilner Scissors and Adson Forceps, single pack 20 units
- SPKA1004 Suture Pack inc Kilner Scissors, Adson Forceps, Mayo Hegar & No3 Handle, double packed 10 units
- SPKA1005 Trauma & Minor Injury Suture Pack, double packed 10 units
- SPKM2003-S Suture Pack inc Kilner Scissors, McIndoe Forceps & Mosquito, single packed 20 units
- FSP2002-S Fine Suture Pack, single packed 20 units.



Forceps Bipolar

STANDARD BIPOLAR FORCEPS

- SU771015044 Standard Bipolar Forceps 1.0mm plain x 150mm, 3M, 4.0plgs 20 units
- SU771520044 Standard Bipolar Forceps 1.5mm plain x 200mm, 3M, 4.0plgs 20 units
- SU772015044 Standard Bipolar Forceps 2.0mm plain x 150mm, 3M, 4.0plgs 20 units
- SU772020044 Standard Bipolar Forceps 2.0mm plain x 200mm, 3M, 4.0plgs 20 units

ADSON BIPOLAR FORCEPS


- SU791015044 Adson Bipolar Forceps 1.0mm plain x 150mm, 3M, 4.0plgs 20 units

BAYONET BIPOLAR FORCEPS

- SU781020044 Bayonet Bipolar Forceps 1.0mm plain x 200mm, 3M, 4.0plgs 20 units
- SU782020044 Bayonet Bipolar Forceps 2.0mm plain x 200mm, 3M, 4.0plgs 20 units

BIPOLAR CABLE

- SU11443075 3M Bipolar 4.0mm Plugs UK/Euro Socket, 20 units
- SU11443076 3M Bipolar 4.0mm Plugs US Socket, 20 units



General & Vascular

● CLAMP COVERS

- Soft silicone jaw covers ensure good grip with reduced trauma
- Supplied in three sizes to suit range of procedures
- Custom sterile packs cut to preferred lengths available on request



STERILE CLAMP COVERS

- CCR2031 Clamp Cover Red (Small) 50mm, double packed 20 units
- CCW2041 Clamp Cover White (Med) 120mm, double packed 20 units
- CCB2051 Clamp Cover Blue (Large) 120mm, double packed 20 units

CLAMP COVERS Supplied non-sterile

- SM3030 Red - Small 5 x 1m
- SM3040 White - Medium 5 x 1m
- SM3050 Blue - Large 5 x 1m

Includes green colour

● SILICONE SLINGS

- Range of pairs and multipacks in 1.5mm and 2.0mm sizes
- New 5mm red & blue for large vessel procedures
- Green colour 2mm available
- Atraumatic oval cross section for less damage during surgery
- Radio opaque and biocompatible medical grade silicone
- Ideal for vessel and nerve retraction, identification, temporary ligation and occlusion



2.0mm Slings

- SM2002 Multicoloured Slings 2.0mm packed in 5's 20 units
- SM2008 Silicone Sling Multipack 2.0mm, (3xRed, 3xBlue, 1xYellow, 1xWhite), double packed 20 units
- SM2030 Red Slings 2.0mm 25 units
- SM2040 Blue Slings 2.0mm 25 units
- SM2050 Yellow Slings 2.0mm 25 units
- SM2060 White Slings 2.0mm 25 units



1.5mm Slings

- SM2100 Multicoloured Slings 1.5mm packed in 5's 20 units
- SM2103 Red Slings 1.5mm 25 units
- SM2104 Blue Slings 1.5mm 25 units
- SM2105 Yellow Slings 1.5mm 25 units
- SM2106 White Slings 1.5mm 25 units

5.0 mm Slings

- SM2513 Silicone Slings Red 5mm, individually double packed 25 units
- SM2514 Silicone Slings Blue 5mm, individually double packed 25 units

TIBBS & OTHER ARTERIAL CANNULA


- Individually packed to avoid wastage or extra reprocessing
- Easier to use than reprocessing very fine tube lumens
- High quality stainless steel suitable for clinicians to bend to shape
- Removes contamination risks through flushing
- Full range to suit procedure, including large vessels

Standard Sizes

- TAC20LC Tibbs Arterial Cannula Large Cone, std length, double packed 20 units
- TAC20SC Tibbs Arterial Cannula Small Cone, std length, double packed 20 units
- TAC20LO Tibbs Arterial Cannula Large Olive, std length, double packed 20 units
- TAC20SO Tibbs Arterial Cannula Small Olive, std length, double packed 20 units

Special Sizes

- TAC20MO Tibbs Arterial Cannula Micro Olive, 2.5mm OD, std length, double packed 20 units
- TAC20XLC Tibbs Arterial Cannula X Large Cone 8mm OD, extra long, double packed 20 units
- TAC20XXLC Tibbs Arterial Cannula XX Large Cone 12.25mm OD, extra long, double packed 20 units
- TAC2003 Tibbs Arterial Cannula, Set of 3 sizes, double packed 20 units
- HFN2001 Heparin Flushing Needle, double packed 20 units




Excellent value

VASCULAR BOOTIES

- Sterile single-use instrument jaw covers, (shods)
- Soft jaw covers provide improved grip for delicate sutures
- Longer length ensure better fit on jaws
- Convenient packed in 6's or 10's
- Alternative colours and sizes available (see Tip Shields for options)
- Range for straight or curved artery forceps and vascular clamps

VB5080 Vascular Booties, single packed in 6's 50 units
 VB5080-D Vascular Booties, double packed in 6's 50 units
 VB5080/10 Vascular Booties, single packed in 10's 50 units
 VB5190 Vascular Booties, double packed in 6's 50 units



Vascular Booties

VASCULAR MOSQUITOS

- Sterile single-use mosquitos with booties ready for use
- Quick and easy to use
- Convenient pack of three ready shod forceps
- Booties fitted correctly everytime

VMQT2003 Vascular Mosquito Forceps with Booties x 3, double packed 20 units




Vascular Mosquito Forceps

Non scalpel procedure

VASECTOMY SET

- High quality stainless steel
- Sterile Single-use avoids cross contamination
- Particularly suitable where local re-processing absent
- Readily available to minimise waiting time

VAS1002 Vasectomy set, double packed 10 units




Vasectomy set

Probes

● VNUS CLOSURE & PHLEBECTOMY HOOKS

VNS1011	VNUS Closure Set with Hook size 1 (2.0mm) , double packed 10 units
VNS1022	VNUS Closure Set with Hook size 2 (1.5mm) , double packed 10 units
VNS1033	VNUS Closure Set with Hook size 3 (1.0mm) , double packed 10 units
PBH2001	Phlebectomy Hook Size 1 (2.0mm), double packed 20 units
PBH2002	Phlebectomy Hook Size 2 (1.5mm), double packed 20 units
PBH2003	Phlebectomy Hook Size 3 (1.0mm), double packed 20 units



NEW

● STOMA BRIDGE

- Three small sizes to suit modern practice
- Easier for patients to wear
- Simple screw thread for quick positioning

STB2040 Stoma Bridge 40mm, double packed, 20 units

STB2050 Stoma Bridge 50mm, double packed, 20 units

STB2060 Stoma Bridge 50mm, double packed, 20 units



Suction Handles

The DTR range of Sterile Single-use Silicone Tubing is recommended for use with these suction handles. For details please refer to the General & Vascular Catalogue

BELLUCCI SUCTION HANDLES

- Metal construction with suction control
- Removes need to clean and re-process the long narrow lumen
- Suitable for theatre use and high risk procedures

BLS2001 Bellucci Suction Handle 175mm x 2mm OD (6fg), double packed 20 units



NEW Range

LONG SUCTION HANDLES

- Metal construction with suction control
- Effective alternative to re-processing very long narrow lumen
- Range of sizes from 280mm x 6fg to 550mm x 12fg
- Suitable for theatre use and high risk procedures

280mm

BLS2280 Micro Laryngeal Suction Handle 280mm x 2mm OD (6fg), double packed 20 units

330mm

BLS2330 Micro Laryngeal Suction Handle 330mm x 2mm OD (6fg), double packed 20 units

BLSF2339 Micro Laryngeal Suction Handle 330mm x 3mm OD (9fg), double packed 20 units

360mm

BLSF2369 Micro Laryngeal Suction Handle, 360mm x 3mm OD (9fg), double packed 20 units



450mm

BLS2450 Bronchoscopy Suction Handle 450mm x 4mm OD (12fg), double packed 20 units


550mm

BLS2550 Oesophageal Suction Handle 550mm x 4mm OD (12fg), double packed 20 units

BLS2550 Sigmoidoscope Suction Handle 560mm x 6mm OD (18fg), double packed 20 units



INSTRUCTIONS FOR USE Frazier Suction Handles and bending for user preference



All dimensions are in mm

1. The recommended points at which to bend the device from each end. Ensure the suction control is in the favoured orientation.



3. Apply firm force by pushing from the thumbs up against the index fingers bending the tube but ensuring not to exceed an angle of 40 degrees



2. Place both thumbs at the point of which requires the bend, and place index fingers firmly gripping the shaft for stability.



4. The top image shows the product bent at the recommended distance from the handle, while the bottom image shows the bend from the tip.

NEW Range

FRAZIER SUCTION HANDLES

- Strong construction suitable for bending to shape during procedure
- Single-use avoids the need to re-straighten or risk breakage
- Four tube diameters available
- Tip options include non-fenestrated, fenestrated with two holes and cross fenestrated with four holes

2.00mm (6fg)

- FSF1001 Frazier Suction Handle 180mm straight fenestrated, double packed 20 units
- FSH2001 Frazier Suction Handle 2mm OD (6fg) 200mm straight, double packed 20 units
- FSF2006 Frazier Suction Handle 2mm OD (6fg) 200mm straight, fenestrated, double packed 20 units
- FSF2046 Frazier Suction Handle 2mm OD (6fg), straight, Cross fenestrated, double packed 20 units

2.66mm (8fg)


- FSF2008 Frazier Suction Handle 2.66mm OD (8fg), straight, fenestrated, double packed 20 units
- FSF2048 Frazier Suction Handle 2.66mm OD (8fg), straight, Cross fenestrated, double packed 20 units
- FSH2008 Frazier Suction Handle 2.66mm OD (8fg), straight, double packed 20 units

3.00mm (9fg)

- FSH2009 Frazier Suction Handle 3mm OD (9fg), 200mm straight, double packed 20 units
- FSF2009 Frazier Suction Handle 3mm OD (9fg), 200mm straight, fenestrated, double packed 20 units
- FSF2049 Frazier Suction Handle 3mm OD (9fg), 200mm straight, Cross fenestrated, double packed 20 units

4.00mm (12fg)

- FSH2012 Frazier Suction Handle 4mm OD (12fg), 200mm double packed 20 units
- FSF2012 Frazier Suction Handle 4mm OD (12fg), 200mm fenestrated, double packed 20 units
- FSF2412 Frazier Suction Handle 4mm OD (12fg), 200mm Cross fenestrated, double packed 20 units



SUCTION REGULATOR

- House style regulator with control
- Luer connector suitable for Verhoeven and other micro suction tubes
- Suitable for theatre use and high risk procedures

DSRP5001 Suction Regulator, black, double packed 50 units



Suction Regulator

MICRO SUCTION TUBES

- Stainless steel tube with colour coded luer connector for micro suction use
- Also suitable for use with syringe as irrigation device
- Particularly suitable for theatre use and high risk procedures
- Available in 86mm, 110mm & 120mm lengths

86mm

MSD2001 Micro Suction Device 86mm x 2mm OD (6fg), double packed 20 units

110mm

MSV50114 Micro Suction (Verhoeven) Green 110mm x 14g (2.10mm OD), double packed 50 units

MSV50117 Micro Suction (Verhoeven) White 110mm x 17g (1.47mm OD), double packed 50 units

MSV50119 Micro Suction (Verhoeven) Brown 110mm x 19g (1.05mm OD), double packed 50 units

MSV50122 Micro Suction (Verhoeven) Black 110mm x 22g (0.70mm OD), double packed 50 units

MSV50127 Micro Suction (Verhoeven) Clear 110mm x 27g (0.40mm OD), double packed 50 units

120mm


MSV50214 Micro Suction (Verhoeven) Green 120mm x 14g (2.10mm OD), double packed 50 units

MSV50217 Micro Suction (Verhoeven) White 120mm x 17g (1.47mm OD), double packed 50 units

MSV50219 Micro Suction (Verhoeven) Brown 120mm x 19g (1.05mm OD), double packed 50 units

MSV50222 Micro Suction (Verhoeven) Black 120mm x 22g (0.70mm OD), double packed 50 units

MSV502127 Micro Suction (Verhoeven) Clear 120mm x 27g (0.40mm OD), double packed 50 units



Micro Suction (Verhoeven)

Complete Range

● POOLE SUCTION HANDLE

- Stainless steel for balance and precision
- Supplied with detachable outer to spread suction
- Ideal when plastic handles are inappropriate
- 240mm in length
- Available in three lengths 190mm 240mm, 330mm

PSHN1001 Poole Suction Handle Neo-Natal, double packed 10 units

PSH1001 Poole Suction Handle, double packed 10 units

PSHL1001 Poole Suction Handle Long, double packed 10 units



Choice of Roses

● YANKAUER SUCTION HANDLES

- Stainless steel for balance and precision
- Ideal when plastic handles are inappropriate
- Detachable rose with three or six holes to spread suction
- Ideal for tonsillectomy
- 305mm length

YST1000 Yankauer Suction Handle 305mm with rose, double packed 10 units

YST1000R Yankauer Suction Handle 305mm with modified rose, double packed 10 units



Tip Shields

- Wide choice of colour size and shape
- Cost effective instrument protection
- Protects staff during re-processed and sterilisation

- **Supplied non sterile**

SM5010	Flat 26mm width, Black 100 units
SM5020	Flat 13mm width, Brown 100 units
SM5030	Flat 6.5mm width, Orange 100 units
SM5040	Round 10mm diameter, Mauve 100 units
SM5050	Round 6mm diameter, Yellow 100 units
SM5060	Round 4mm diameter, Red 100 units
SM5070	Round 3mm diameter, Green 100 units
SM5080	Round 2mm diameter, Blue 100 units
SM5090	Round 1mm diameter, White 100 units



Tubing Products

- **SILICONE TUBING**

- Full range of wall and lumen sizes
- Wide application of uses
- Supplied in 15m packs for cutting to size
- Suitable for custom procedure trays

- **Supplied non sterile**

SM105	ID 0.3	OD 0.64	15m roll
SM135	ID 0.51	OD 0.94	15m roll
SM155	ID 0.64	OD 1.19	15m roll
SM175	ID 0.76	OD 1.65	15m roll
SM205	ID 1.02	OD 2.16	15m roll
SM235	ID 1.47	OD 1.96	15m roll
SM265	ID 1.57	OD 2.41	15m roll
SM285	ID 1.57	OD 3.18	15m roll
SM305	ID 1.98	OD 3.18	15m roll
SM305P	ID 1.98	OD 3.18	15m roll Blue Pigmented
SM325	ID 2.65	OD 4.88	15m roll
SM335	ID 3.35	OD 4.65	15m roll
SM365	ID 3.18	OD 6.35	15m roll
SM405	ID 4.80	OD 7.90	15m roll
SM415	ID 5.00	OD 8.30	15m roll
SM425	ID 4.80	OD 9.50	15m roll
SM445	ID 6.40	OD 9.50	15m roll
SM465	ID 6.40	OD 11.10	15m roll
SM485	ID 6.40	OD 12.70	15m roll
SM505	ID 7.90	OD 12.70	15m roll
SM525	ID 9.50	OD 12.70	15m roll
SM545	ID 9.50	OD 14.30	15m roll
SM565	ID 9.50	OD 15.90	15m roll
SM625	ID 12.70	OD 17.50	15m roll
SM645	ID 12.70	OD 19.05	2x7.5m roll
SM685	ID 15.90	OD 22.23	2x7.5m roll



New Extended Range

STERILE SILICONE TUBE SETS

- Sterile for extra confidence and where lengths over 270mm are impossible to validate
- Cost effective alternative to in-house production
- Full range of tube sizes cut to consistent lengths
- Complimentary to suction handle range
- Tailored to suit needs including mixed packs

DST105

DST10503 Sterile Silicone Tube 0.3 ID x 0.64mm OD x 300mm long, double packed 10 units

DST175

DST17505 Sterile Silicone Tube 0.76 ID x 1.65mm OD x 515mm long, double packed 10 units

DST205

DST2051004 Sterile Silicone DCR Tube 1.02mm ID x 2.16mm OD x 40mm long, double packed 10 units

DST20503 Sterile Silicone Tube 1.02mm ID x 2.16mm OD x 300mm long, double packed 10 units

DST20510 Sterile Silicone Tube 1.02mm ID x 2.16mm OD x 1000mm long, double packed 10 units

DST285

DST2851 Sterile Silicone Tube 1.57mm ID x 3.1mm OD x 250mm long, double packed 10 units

DST2852 Sterile Silicone Tube 1.57mm ID x 3.1mm OD x 250mm long two pieces, double packed 10 units

DST28503 Sterile Silicone Tube 1.57mm ID x 3.18mm OD x 300mm, long double packed 10 units

DST28504 Sterile Silicone Tube 1.57mm ID x 3.1mm OD x 450mm long, double packed 10 units

DST28507 Sterile Silicone Tube 1.57mm ID x 3.1mm OD x 700mm long, double packed 10 units

DST28530 Sterile Silicone Tube 1.57mm ID x 3.18mm OD, 3000mm long, double packed 20 units

DST305

DST305002 Sterile Silicone Tube 1.98mm ID x 3.18mm OD x 25mm long, double packed 10 units

DST30502 Sterile Silicone Tube 1.98mm ID x 3.18mm OD x 200mm long, double packed 10 units

DST30503 Sterile Silicone Tube 1.98mm ID x 3.18mm OD x 300mm long, double packed 10 units

DST30504P Sterile Silicone Tube 1.98mm ID x 3.18mm OD x 450mm long, blue double packed 20 units

DST325

DST325002 Sterile Silicone Tube 2.65mm ID x 4.88mm OD x 25mm long, double packed 10 units

DST32503 Sterile Silicone Tube 2.65mm ID x 4.88mm OD x 300mm long, double packed 10 units

DST335

DST335815 Sterile Silicone tube 3.35mm ID x 4.65mm OD x 150mm long x 8 pieces double packed 20 units

DST365

DST36502 Sterile Silicone Tube 3.2mm ID x 6.4mm OD x 200mm long, double packed 20 units



- DSTS36502 Sterile Silicone Tube 3.2mm ID x 6.4mm OD x 200mm long with spigot, double packed 20 units
- DST36511 Sterile Silicone Tube 3.2mm ID x 6.4mm OD, 1100mm long, double packed - House 10 units
- DST36508 Sterile Single-use Silicone Tube 3.2mm ID x 6.4mm OD x 800mm long, double packed 10 units
- DST36536 Sterile single-use silicone tube 3.2mm ID x 6.4mm OD x 3600mm long, double packed 10 units

DST405

- DST405015 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 150mm long, double packed 10 units
- DST40503 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 300mm long, double packed 10 units
- DST40506 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 600mm long, double packed 10 units
- DST40576-S Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 760mm long, single packed 10 units
- DST40510 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 1000mm long, double packed 10 units
- DST40512 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 1200mm long, double packed 10 units
- DST40513 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 1350mm long, double packed 20 units
- DST40515 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 1500mm long, double packed 20 units
- DST40517 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 1750mm long, double packed 20 units
- DST4052M Sterile Silicone Tube 4.8mm ID x 7.9mm OD lengths 1 x 1.8m & 1 x 0.45m, double mixed pack 20 units
- DST40525 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 2500mm long, double packed 20 units
- DST40536 Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 3600mm long, double packed 20 units

DST415

- DST41506 Sterile Silicone Tube 5.0mm ID x 8.3mm OD x 600mm long, double packed
- DST41512 Sterile Silicone Tube 5.0mm ID x 8.3mm OD x 1200mm long, double packed
- DST41515 Sterile Silicone Tube 5.0mm ID x 8.3mm OD x 1500mm long, double packed 20 units
- DST41525 Sterile Silicone Tube 5.0mm ID x 8.3 mm OD x 2500mm, double packed 20 units
- DST41536 Sterile Silicone Tube 5.0mm ID x 8.3 mm OD x 3600mm long, double packed 20 units

DST425

- DST42512 Sterile Silicone Tube 4.8mm ID x 9.5mm OD x 1200mm long, double packed
- DSTS42576 Sterile Silicone Tube 4.8mm ID x 9.5mm OD x 760mm long with Spigot, double packed 10 units
- DST4252M Sterile Silicone Tube 4.8mm ID x 9.5mm OD lengths 1 x 1.8m & 1 x 0.45m, double mixed pack 20 units

DST445

- DST44503 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 300mm long, double packed 10 units
- DST44507 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 700mm long, double packed 10 units
- DST44510 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 1000mm long, double packed 10 units
- DSTS44510 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 1000mm long with Spigot, double packed 10 units
- DST44514 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 1450mm long, double packed 10 units
- DST44515 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 1500mm long, double packed 20 units
- DST44525 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 2500mm long, double packed 20 units
- DST44536 Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 3600mm long, double packed 20 units

DST465

- DST46504 Sterile Silicone Tube 6.4mm ID x 11.1mm OD x 450mm long, double packed 20 units
- DST46507 Sterile Silicone Tube 6.4mm ID x 11.1mm OD x 750mm long, double packed 20 units
- DST46510 Sterile Silicone Tube 6.4mm ID x 11.1mm OD x 1000mm long, double packed 20 units

DST485

- DST48525 Sterile Silicone Tube 6.4mm ID x 12.7mm OD x 2500mm long, double packed 10 units

DST505

- DST50515 Sterile Silicone Tube 7.9mm ID x 12.7mm OD x 1500mm long, double packed 10 units

DST525

- DST52536 Sterile silicone tube 9.5mm ID x 12.7mm OD x 3600mm long, double packed 10 units

DST545

- DST545007 Sterile Silicone Tube 9.5mm ID x 14.3mm OD x 75mm long, double packed 10 units
- DSTS54502 Sterile Silicone Tube 9.5mm ID x 14.3mm OD x 250mm long with Spigot, double packed 10 units

DST565

- DST56515 Sterile Silicone tube 9.50mm ID x 15.90mm OD x 1500mm, double packed 10 units

DST625

- DST6250015 Sterile Silicone Tube 12.7mm ID x 17.5mm OD x 15mm long, double packed
- DST62517 Sterile Silicone Tube 12.7mm ID x 17.50mm OD x 1750 long, double packed 10 units
- DST62525 Sterile Silicone Tube 12.7mm ID x 17.50mm OD x 2500 long, double packed 20 units

SILICONE BANDS

- Sterile for extra confidence
- Cost effective alternative to in house production
- Latex free to replace elastic bands
- Can be tailored to suit needs

DSB3211-S Silicone Bands 2.16mmx110mm dia x 3,
single packed 10 units


TUBE ADAPTORS

DBTL1001-S Blue Tube Adapter Long (Catheter Mount
Spigot), single packed 10 units

DBTS1001-S Blue Tube Adapter Short (Catheter Mount
Spigot), single packed 10 units

BLUE CLAMPS

DBC5051-S Blue Clamp Mediplast, single packed
10 x 5 units



Notes



DTR Medical

Our mission is to provide medical professionals with high quality single-use sterile disposables that deliver patient and clinical value combined with exceptional levels of service.

We operate our own production facilities equipped with ISO class 7 cleanrooms and the company is accredited to ISO9001:2008, ISO13485:2003 and MDD2007/47/EC.



INVESTOR IN PEOPLE

DTR Medical Ltd

17 Clarion Court, Enterprise Park, Swansea, SA6 8RF, UK
Telephone: +44 (0)1792 797910 Fax: +44 (0)1792 797955
Email: customerservice@dtrmedical.com

www.dtrmedical.com