


## SCHUKNECHT SUCTION HANDLE **NEW**

Removes the difficulty of cleaning a very narrow lumen

Suitable for theatre use and high risk procedures

Range of sizes from 20 to 24 standard wire gauge



Code	Description	Unit of Sale
SHK2020	Sterile Single-use Schuknecht Suction Handle 20swg, double packed	20 units
SHK2022	Sterile Single-use Schuknecht Suction Handle 22swg, double packed	20 units
SHK2024	Sterile Single-use Schuknecht Suction Handle 24swg, double packed	20 units

## METAL "GLASS-TYPE" SUCTION HANDLE

Metal version of traditional glass suction handle

Suitable for where wide bore is required

Single-use to prevent cross contamination



Code	Description	Unit of Sale
GSH2001	Sterile Single-use Metal "Glass-Type" Suction Handle, 200mmx5mm OD (15fg), double packed	20 units

## HOUSE SUCTION HANDLES

Dual purpose suction and irrigation

Convenient suction control

Complete range of sizes from 0.7mm to 5.0mm

Single-use overcomes repairs to tips damaged through drilling

Particularly suitable for collection of bone dust



Code	Description	Unit of Sale
HST0109/07	Sterile Single-use House Suction & Irrigation Handle (0.9 x 0.7mm), double packed	20 units
HST0112/09	Sterile Single-use House Suction & Irrigation Handle (1.2 x 0.9mm), double packed	20 units
HST0115/12	Sterile Single-use House Suction & Irrigation Handle (1.5 x 1.2mm), double packed	20 units
HST0118/18	Sterile Single-use House Suction & Irrigation Handle (1.8 x 1.8mm), double packed	20 units
HST0125/20	Sterile Single-use House Suction & Irrigation Handle (2.5 x 2.0mm), double packed	20 units
HST0130/20	Sterile Single-use House Suction & Irrigation Handle (3.0 x 2.0mm), double packed	20 units
HST0145/20	Sterile Single-use House Suction & Irrigation Handle (4.5 x 2.0mm), double packed	20 units
HST0145/30	Sterile Single-use House Suction & Irrigation Handle (4.5 x 3.0mm), double packed	20 units
HST0145/45	Sterile Single-use House Suction & Irrigation Handle (4.5 x 4.5mm), double packed	20 units
HST0150/45	Sterile Single-use House Suction & Irrigation Handle (4.5 x 4.5mm), double packed	20 units

## FENESTRATED HOUSE SUCTION HANDLE

Dual purpose suction and irrigation

Convenient suction control

Fenestrated tip to spread suction

Single-use overcomes repairs to tips damaged through drilling



Code	Description	Unit of Sale
HSLF118/18	Sterile Single-use House Fenestrated Handle (1.8 x 1.8mm), double packed	20 units